High purity laboratory consumables for advanced sciences.



EXAPURE TM

Certified extractable and metal-free.

High purity laboratory consumables for advanced sciences.



DISPOSABLE FILTRATION DEVICES

Syringe filters, disc membrane filters, in-line filters

Whenever there is a need to filter a viscous solution, a high particulate medium or simply a precious sample, EXAPURÉ™ Filters provide easy, effective, and high purity devices.

> No need to hesitate: EXAPURE™.

SPECIAL FEATURES

- · High purity / HPLC certified High purity polymers (USP grade) assure that the filter will not create artifacts. Glues and sealants are excluded to minimise extractables.
- Wide chemical compatibility High chemical resistance for broader measurement methods.
- Optimised design minimises dead volume and optimised membrane structure reduces blockage and increases volume throughput.
- Filtration efficiency The microglass prefilter reduces blockage and increases flow.

· Filtration efficiency & design · Leak tight / Proven integrity

All connectors match with a variety of syringes and tubings. Each lot is 100% tested to ensure product integrity and reliability.

Proven compatibility Devices are manufactured to standard dimensions.



13 mm diameter



25 mm diameter



sterile, 25 mm diameter individually packed

SYRINGE FILTERS The EXAPURE™ Syringe Filters are specifically designed for purification, isolation and separation. The disposable syringe filters are made of carefully selected membranes to ensure the highest filtration reliability. Moreover, the filter design is optimised for the highest flow while minimising the dead volume. They are made with of a wide variety of different membranes with a polypropylene housing using the most advanced methods and design features available today. They are particularly appreciated for various applications in pharmaceutical, environmental, biotechnology, food, beverage and water analysis laboratories, among others.

> The EXAPURE™ Syringe Filters are exclusively made of a high purity polymer housing and heat sealed without the use of glues or any ensuring absolute retention and extremely low The housing is made of a specifically selected polypropylene grade for sterile applications as well as against highly corrosive chemicals. All EXAPURE™ Syringe Filters are HPLC certified and metal-free.

The EXAPURE™ Syringe Filters PA are made of a high chemical resistant nylon (polyamide) membrane, which has been designed for applications as sensitive as the filtration of HPLC solutions.

The EXAPURE™ Syringe Filters PTFE are made of a high chemical resistant Teflon (polytetrafluoroethylene) membrane, which has been designed for applications as sensitive as the filtration of HPLC and GC solutions.

The EXAPURE™ Syringe Filters PES have been designed for the filtration of sterile media, serums, buffers and other aqueous solutions.
The glass microfiber prefilter ensures an economical filtration of serums while maintaining the filtration capacity of the PES membrane, which is specifically selected for its low protein binding properties. In addition, the individual packaging ensures the filter sterility until utilisation. EXAPURE™ Syringe Filters PES are pyrogen and metal-free.

Technical Specifications

		40	0.5		
Membrane diameter (mm)	4	13	25		
Effective filtration area (cm²)	0.07	0.8	4.15		
Hold-Up volume with purge (µI)	< 5	< 20	< 50		
Max. operating pressure	≥ 5 kg/cm² @ 30°C				
Biosafety	All materials pass USP Class VI test for plastics				
Materials of construction	Virgin Polypropylene housing with pure and metal- free membrane				
Connectors	Female Luer Lock (FLL) Male Slip Luer (MSL)				
Flow direction	Flow from inlet to outlet (FLL to MSL)				
Sterilisation	ETO or autoclavable at 121°C for 30 minutes				

EXAPURE™ vs. Competition						
5	Competitor S					
4 Competitor						
3						
2- EXAPURE™						
1 PES						
Protein binding test (0.75%	BSA, μg/cm ²)					

Article	Reference	Diam./pore	Packaging	Quantity C	olour Ster	ile Note
EXAPURE™ Syringe Filters PA 0.45-13, HPLC	ATSY25NN4	(mm) / (µm) 13 / 0.45)	100 gr	reen no	Ø 13 mm polyamide membrane (Nylon), pore size
EXAPURE™ Syringe Filters PA 0.2-25, HPLC	ATSY25NN2	25 / 0.2	bag	100 bl	lue no	 0.45 μm, non sterile 25 mm polyamide membrane (Nylon), pore size 0.2 μm, non sterile
EXAPURE™ Syringe Filters PA 0.45-25, HPLC	ATSY25NN4	25 / 0.45	bag	100 gr	reen no	 25 mm polyamide membrane (Nylon), pore size 0.45 μm, non sterile
EXAPURE™ Syringe Filters PTFE 0.2-13, GC	ATSY13TF2	13 / 0.2	bag	100 re	ed no	Ø 13 mm PTFE membrane (Teflon), pore size 0.2 um, non sterile, ultra-low extractables
EXAPURE™ Syringe Filters PTFE 0.45-13, GC	ATSY13TF4	13 / 0.45	bag	100 bl	lack no	Ø 13 mm PTFE membrane (Teflon), pore size 0.45 μm, non sterile, ultra-low extractables
EXAPURE™ Syringe Filters PTFE 0.2-25, GC	ATSY25TF2	25 / 0.2	bag	100 re	ed no	Ø 25 mm PTFE membrane (Teflon), pore size 0.2 μm, non sterile, ultra-low extractables
EXAPURE™ Syringe Filters PTFE 0.45-25, GC	ATSY25TF4	25 / 0.45	bag	100 bl	lack no	Ø 25 mm PTFE membrane (Teflon), pore size 0.45 μm, non sterile, ultra-low extractables
EXAPURE™ Syringe Filters PTFE 1.0-25, GC			bag	,	yan no	Ø 25 mm PTFE membrane (Teflon), pore size 1.0 μm, non sterile, ultra-low extractables
EXAPURE™ Syringe Filters PES 0.2-25, PF-GF, sterile	ATSY25PGS2	25 / 0.2	individual	100 w	hite yes	Ø 25 mm polyethersulfone membrane (PES), por size 0.2 μm, ETO sterilisation, with glass microfiber prefilter, pyrogen and metal-free
EXAPURE™ Syringe Filters PES 0.45-25, PF-GF, sterile	ATSY25PGS4	25 / 0.45	individual	100 w	hite yes	Ø 25 mm polyethersulfone membrane (PES), por size 0.45 µm, ETO sterilisation, with glass microfiber prefilter, pyrogen and metal-free
EXAPURE™ Syringe Filters GF 1.0-25, GC	ATSY25GF10	25 / 1.0	bag	100 w	hite no	Ø 25 mm glass microfiber membrane, pore size 1.0 μm, non sterile, ultra-low extractables



47 mm diameter



sterile, 47 mm diamete individually packed

MEMBRANE DISC FILTERS The EXAPURE™ Disc Filters bring to the laboratory user a range of membrane filters whose advanced technical specifications ensure today's most advanced filtration performance for a wide range of applications. The membrane disc filters offer accurately controlled pore size distribution and higher strength and flexibility which ensure reproducibility and consistency. Moreover, the membrane material is carefully selected to release the lowest extractable. Sterile and individually packed membranes disc filters as well as coloured and gridded types are available for specialised applications. EXAPURE™ Disc Filters pass the biological test for Class VI plastics as described in USP and the oxidizable matters are systematically controlled to be within limits specified in USP. All EXAPURE™ Disc Filters are metal-free.

Technical Specifications

recillical opecifications						
Pore size (µm)	0.2 0.45 0.8					
Water flow rate 10 psi / 20°C (ml/min.cm²)	≥ 14.5	≥ 38.5	≥ 125			
Retention efficiency (LRV/cm²)	≥ 7	≥ 7	≥ 7			
	for Br. diminuta	for Sr. marcescens	for yeast			
Max. operating pressure	5 kg/cm ² @ 30°C					
Biosafety	All materials pass USP Class VI test for plastics					
Materials of construction	High purity and metal-free membrane					
Sterilisation	ETO or autoclavable at 121°C for 30 minutes					

Reference Diam./pore Packaging Quantity Colour Sterile Not (mm) / (µm) 47 / 0.2 box EXAPURE™ Disc Filters PA 0.2-47 AT47HNN2 Ø 47 mm polyamide membrane (Nylon), pore size 0.2 μm Ø 47 mm polyamide membrane (Nylon), pore size 0.45 μm Ø 47 mm cellulose nitrate membrane, pore size EXAPURE™ Disc Filters PA 0.45-47 47 / 0.45 box Ø 47 mm cellulose nitrate membrane, pore size 0.45 µm, ETO sterilisation, gridded EXAPURE™ Disc Filters CN 0.45-47, gridded. 47 / 0.45 individual sterile EXAPURE™ Disc Filters Pack CN212, 100 ml Funnels, CN 0.45-47, gridded, sterile ATFM47GCNS4 47 / 0.45 individual 212 x \varnothing 47 mm cellulose nitrate membrane, pore size 0.45 μ m, ETO sterilisation, gridded 212 x absorbent pad 12 x autoclavable 100 ml PP funnel with cover



50 mm diameter

IN-LINE FILTERS The EXAPURE™ In-Line Filters feature a housing made out of a specifically selected high purity polypropylene grade to maintain sample purity and for sterile applications as well as against highly corrosive chemicals. They are available with a choice of filtration membrane to suit a broad range of applications with aqueous, gaseous and organic samples. The most advanced technologies have been used to ensure convenience and performance for inline sterilisation of gases and bacterial air venting in various pharmaceutical and biotechnology processes. The large 50 mm in-line filters are designed for larger volume sample filtration in the laboratory, at a pilot plant or in manufacturing. The EXAPURE™ In-Line Filters pass the biological test for Class VI plastics as described in USP, the extractables with isopropyl alcohol and oxidizable matters are both systematically controlled to be within limits specified in USP. IPF Filters are specifically designed to connect directly into an HPLC line to simultaneously filter and degas the mobile phase as it is being used. Smaller in-line filters in 37 mm and 25 mm diameters are also available for specific applications. ITF Filters are outstanding for sterilisation or air venting applications, are microbially validated as per ASTM 838-83, complies with USFDA 21 CFR 211.72 and USDFA 21 CFR 177.1520. All EXAPURE™ In-Line Filters are autoclavable, individually packed and metal-free.

Technical Specifications

recinical opecinications						
Membrane diameter (mm)	25	37	50			
Effective filtration area (cm ²)	4.15	9.6	18			
Max. operating pressure	≥ 3 kg/cm ² @ 30°C					
Biosafety	All materials pass USP Class VI test for plastics					
Materials of construction	Virgin Polypropylene housing with pure and metal- free membrane					
Connectors	Stepped hose barb					
	Female Luer Lock (FLL) Male Slip Luer (MSL)					
Flow direction	Flow from inlet to outlet (FLL to MSL)					
Sterilisation	ETO or autoclavable at 121°C for 20 minutes					

The EXAPURE™ In-Line Filters ITF feature a PTFE (Teflon) membrane which is specifically designed for chemically aggressive solutions, reagents and organic solvents. This lightweight device is particularly suitable for protective vents and for in-line filtration and

isolation purposes.

The EXAPURE™ Disc Filters PA are made of

The EXAPURE™ Disc Filters PA are made of a high chemical resistant nylon (polyamide) membrane, which has been designed for applications as sensitive as the filtration of HPLC solutions. They are are hydrophilic, non-media migrating, biologically inert, designed for aqueous as well as organic solvent filtration. EXAPURE™ Disc Filters PA

specifically developed for water microbiology and validated to ASTM D4200-82 and D3508-

78. They are made of a high grade cellulose nitrate membrane, hydrophilic, presterilised

and individually packed in a form of roll which is compatible with most disc membrane

■ EXAPURE™ Disc Filters PA

■ EXAPURE™ Disc Filters CN


0.2 0.45 0.8 1.2 Water flow rate (ml/min·cm²) vs. pore size (μm)

are HPLC certified and metal-free. The EXAPURE™ Disc Filters CN gridded are

dispenser

300-200

The EXAPURE™ In-Line Filters IPF feature a PP (polypropylene) which provides excellent flow rate characteristics and is ideal for filtering large volumes to 1 litre or more of aqueous and solvent samples. They possess a high-flow polypropylene membrane for mobile phases containing organic solvents or other highly corrosive liquids such as acids or bases The 1/8" MNPT filter connection is designed to match the most difficult specifications for leak-free applications. The EXAPURE™ In-Line Filters IPF filter with ease highly loaded as well as viscous solutions without leaking, unlike filters with horse-barb connections. The 5 μm pore size membrane filters particularly well the coarser particles that would rapidly obstruct conventional 0.2 or 0.45 µm filters. The EXAPURE™ In-Line Filters PP 5-50 fulfil the most demanding requirements in the chemical, but also in the food and ink industry.



Air flow rate (dm3/min) vs. diameter (mm)

ALYS Technologies Labware business unit

e**chnologies SA** Bas-de-Plan 11 CH-1030 Bussigny-près-Lausanne Tel : +41 (21) 312 42 60 Fax : +41 (21) 312 42 61 labware@alys-technologies.com www.alys-technologies.com

Distributed by:

Article	Reference		Packaging	Quantit	y Colour	Sterile	Note
		(mm) / (μm					
EXAPURE™ In-Line Filters ITF 0.2-50, sterile	AT50ITF2	50 / 0.2	bag	10	white	yes	Ø 50 mm PTFE membrane (Teflon), pore size 0.2
							μm, hose barb connection, ETO sterilisation
EXAPURE™ In-Line Filters ITF 0.4-50, sterile	AT50ITF4	50 / 0.45	bag	10	white	yes	Ø 50 mm PTFE membrane (Teflon), pore size
							0.45 μm, hose barb connection, ETO sterilisation
EXAPURE™ In-Line Filters IPF 5-50	AT50IPF50	50 / 5	bag	10	white	no	Ø 50 mm Polypropylene membrane (PP), pore
							size 5 µm, 1/8" MNPT connection

EXAPURE™ High purity laboratory consumables for advanced sciences.



Reaching your science. Your choice.

EXAPURE™

High purity laboratory consumables info@exapure.com www.exapure.com

EXAPURE™, a Swiss brand

Pure Heart. Pure Mind. Pure Science.